

## Calf implants CE 0483


### POLYsmooth™

Filled with silicone gel, individually packed, sterile

Order N°	Dimensions (cm)			Volume (ml)
	A	B	C	


#### Symmetrical design (Moderate profile)

10200-012 <sup>1</sup>	2.1	6.2	1.6	12
10200-030 <sup>1</sup>	3.0	10.2	1.9	30
10200-070	3.3	15.7	2.3	70
10200-090	3.2	19.5	2.5	90
10200-120	3.4	26.8	2.5	120
10200-140 <sup>1</sup>	4.4	26.2	2.5	140



#### Anatomical design (Anatomical profile)

10201-085	4.7	15.1	2.4	85
10201-140	5.4	19.7	2.5	140
10201-180 <sup>1</sup>	5.8	23.1	2.6	180




## Pectoralis implants CE 0483

Filled with silicone gel, individually packed, sterile

Order N°	Dimensions (cm)			Volume (ml)
	A	B	C	

### POLYsmooth™

10500-190L <sup>1</sup>	14.7	10.2	2.4	190
10500-190R <sup>1</sup>	14.7	10.2	2.4	190
10500-230L <sup>1</sup>	16.0	10.7	2.7	230
10500-230R <sup>1</sup>	16.0	10.7	2.7	230
10500-300L <sup>1</sup>	16.9	11.6	2.9	300
10500-300R <sup>1</sup>	16.9	11.6	2.9	300



### MESMO®

15500-190L <sup>1</sup>	14.7	10.2	2.4	190
15500-190R <sup>1</sup>	14.7	10.2	2.4	190
15500-230L <sup>1</sup>	16.0	10.7	2.7	230
15500-230R <sup>1</sup>	16.0	10.7	2.7	230
15500-300L <sup>1</sup>	16.9	11.6	2.9	300
15500-300R <sup>1</sup>	16.9	11.6	2.9	300

L = Left / R = Right

## Scar treatment - POLYTECH Silicone gel & sheets CE

For the treatment of hypertrophic scars, keloids and erythema. Individually packed, non sterile

Order N°	Specification	Thickness (mm)	Length x Width (mm)
----------	---------------	----------------	---------------------

### Silicone gel

1300-015	Tube, 15 ml		
----------	-------------	--	--

### Silicone gel sheets

1310-010	Square	0.9	100 x 100
1320-005	Rectangle	0.9	200 x 50
1330-002	Rectangle	0.9	300 x 20
1330-025	Rectangle	0.9	300 x 250
1370-070	① Disc (areola)	0.9	70 x 70
1370-071 <sup>1</sup>	② Keyhole (reduction)	0.9	70 x 70 + 35 x 90
1395-340 <sup>1</sup>	③ Anchor (reduction)	0.9	95 x 340

